


Product Datasheet

Mouse BID ELISA Kit (orb865341)

Description	Mouse BID ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-D195
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
Uniprot ID	P70444
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse BID
PicoKine
ELISA Kit
Standard
Cu...