



Product Datasheet

CD11c/Itgax Antibody (orb865421)

Description	CD11c/Itgax Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E.coli-derived mouse CD11c/Itgax recombinant protein (Position: P140-T342).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	127 kDa
Uniprot ID	Q9QXH4
Expiration Date	12 months from date of receipt.



Western blot analysis of CD11c/Itgax using anti-CD11c/Itgax antibody.



IHC analysis of CD11c/Itgax using anti-CD11c/Itgax antibody.



IHC analysis of CD11c/Itgax using anti-CD11c/Itgax antibody.