



Product Datasheet

GPR54/KISS1R Antibody (orb865425)


Description	GPR54/KISS1R Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, WB
Immunogen	E.coli-derived human GPR54/KISS1R recombinant protein (Position: N60-A373).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	80 kDa
Uniprot ID	Q969F8
Expiration Date	12 months from date of receipt.



Western blot analysis of GPR54/KISS1R us...



IHC analysis of GPR54/KISS1R using anti-...



IHC analysis of GPR54/KISS1R using anti-...