



## Product Datasheet

### RNF20 Antibody (orb865457)


|                            |                                                                                                                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | RNF20 Antibody                                                                                                                                                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                 |
| <b>Reactivity</b>          | Human, Rat                                                                                                                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                           |
| <b>Tested Applications</b> | ELISA, FC, ICC, IF, WB                                                                                                                                                                                                                 |
| <b>Immunogen</b>           | E.coli-derived human RNF20 recombinant protein (Position: E123-S389).                                                                                                                                                                  |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                            |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                                                              |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                  |
| <b>Application notes</b>   | Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml. |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                             |
| <b>MW</b>                  | 114 kDa                                                                                                                                                                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q5VTR2</a>                                                                                                                                                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                        |



Western blot analysis of RNF20 using ant...



IF analysis of RNF20 using anti-RNF20 an...



Flow Cytometry analysis of U87 cells usi...