



## Product Datasheet

### Golgin 97/GOLGA1 Antibody (orb865525)


<b>Description</b>	Golgin 97/GOLGA1 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, ICC, IF, IHC, WB
<b>Immunogen</b>	E.coli-derived human Golgin 97/GOLGA1 recombinant protein (Position: M1-K752).
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	88-97 kDa
<b>Uniprot ID</b>	<a href="#">Q92805</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Golgin 97/GOLGA1...



IHC analysis of Golgin 97/GOLGA1 using a...



IHC analysis of Golgin 97/GOLGA1 using a...