
Product Datasheet

**Prothrombin Antibody (monoclonal, 4G13G9)
(orb865671)**

Description

Prothrombin Antibody (monoclonal, 4G13G9)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

E. coli-derived human Prothrombin recombinant protein (Position: Y97-R124).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality


Monoclonal

MW


70 kDa

Uniprot ID
[P00734](#)
Expiration Date


12 months from date of receipt.



Western blot analysis of Prothrombin using anti-Prothrombin antibody.



IHC analysis of Prothrombin using anti-Prothrombin antibody.



IHC analysis of Prothrombin using anti-Prothrombin antibody.