


## Product Datasheet

### His-Tag Mouse Monoclonal Antibody (orb865747)

<b>Description</b>	His-Tag Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, IP, WB
<b>Immunogen</b>	A synthesized peptide
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	WB 1:1000-1:5000 ICC/IF 1:100-1:200 IP 1:100
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of His-tagged prot...