


Product Datasheet

HA-Tag Mouse Monoclonal Antibody (orb865749)

Description	HA-Tag Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthesized peptide
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:1000-1:5000 ICC/IF 1:100-1:200 IP 1:100
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of HA-tagged prote...