


Product Datasheet

GFP-Tag Mouse Monoclonal Antibody (orb865750)

Description	GFP-Tag Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Full-length GFP protein
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:3000-1:10000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:100
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P42212
Expiration Date	12 months from date of receipt.



Western blot analysis of GFP protein exp...