


Product Datasheet

EHMT2/G9A Rabbit Monoclonal Antibody (orb866287)

Description	EHMT2/G9A Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human EHMT2/G9A
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200
Isotype	IgG
Clonality	Monoclonal
MW	160-180 kDa
Uniprot ID	Q96KQ7
Expiration Date	12 months from date of receipt.



Western blot analysis of EHMT2/G9A expression.