


Product Datasheet

TCF21 Rabbit Monoclonal Antibody (orb866324)

Description	TCF21 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	A synthesized peptide derived from human TCF21
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IP 1:60FC 1:80
Isotype	IgG
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	043680
Expiration Date	12 months from date of receipt.



Western blot analysis of TCF21 expression...