


Product Datasheet

Histone H3 (acetyl K27) Rabbit Monoclonal Antibody (orb866691)

Description	Histone H3 (acetyl K27) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from Histone H3 (acetyl K27)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200ChIP 1:50
Isotype	IgG
Clonality	Monoclonal
MW	17 kDa
Uniprot ID	Q6NXT2 , P68431 , P84243 , Q71DI3 , Q16695
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H3 (ace...