

## Product Datasheet

# SARS-CoV-2 SPIKE PROTEIN S2 antibody (orb867039)

<b>Description</b>	Rabbit polyclonal antibody to SARS-CoV-2 SPIKE PROTEIN S2
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Immunogen</b>	Peptide corresponding to 16 amino acids near the center terminus of SARS-CoV-2 (COVID-19) Spike glycoprotein. The immunogen is located within 900-950 amino acids of SARS-CoV-2 (COVID-19) Spike protein.
<b>Target</b>	SARS-CoV-2 SPIKE PROTEIN S2
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Expiration Date</b>	12 months from date of receipt.