

---

## Product Datasheet

**SARS-CoV-2 SPIKE PROTEIN CLEAVAGE SITE  
antibody (orb867042)**

**Description**

Rabbit polyclonal antibody to SARS-

SITE

**Species/Host**

Rabbit

**Reactivity**

Virus

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IF, WB

**Immunogen**

Peptide corresponding to 12 amino acids near the center of SARS-CoV-2 (COVID-19) Spike glycoprotein. The immunogen is located within 650-700 amino acids of SARS-CoV-2 (COVID-19) Spike protein.

**Target**

SARS-CoV-2 SPIKE PROTEIN  
CLEAVAGE SITE

**Preservatives**

Purified IgG - liquid

**Form/Appearance**

Purified IgG - liquid

**Concentration**

1.0 mg/ml

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Purified

**Uniprot ID****PODTC2****Expiration Date**

12 months from date of receipt.