

Product Datasheet

SARS-CoV-2 SPIKE PROTEIN antibody (orb867043)

Description	Rabbit polyclonal antibody to SARS-CoV-2 SPIKE PROTEIN
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Peptide corresponding to 20 amino acids near the carboxy terminus of SARS-CoV-2 (COVID-19) Spike glycoprotein. The immunogen is located within the last 50 amino acids of SARS-CoV-2 (COVID-19) Spike protein.
Target	SARS-CoV-2 SPIKE PROTEIN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.