

Product Datasheet

SARS-CoV NUCLEOPROTEIN antibody (orb867045)

Description	Mouse monoclonal antibody to SARS-CoV NUCLEOPROTEIN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant fragment aa1-49 of the SARS nucleoprotein
Target	SARS-CoV NUCLEOPROTEIN
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P59595
Expiration Date	12 months from date of receipt.