

## Product Datasheet

### CD62P antibody (Biotin) (orb867148)

<b>Description</b>	Mouse monoclonal antibody against CD62P conjugated to Biotin
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	HEK293 derived recombinant human full length CD62P (Met1- Ala771)
<b>Target</b>	CD62P
<b>Preservatives</b>	Purified IgG conjugated to Biotin - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Biotin - liquid
<b>Concentration</b>	IgG concentration 0.1 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P16109</b>
<b>Expiration Date</b>	12 months from date of receipt.