

Product Datasheet

CD258 antibody (Biotin) (orb867160)

Description	Mouse monoclonal antibody against CD258 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	E. coli derived recombinant human CD258 (Leu61-Val240)
Target	CD258
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	O43557
Expiration Date	12 months from date of receipt.