

Product Datasheet

CD130 antibody (Biotin) (orb867164)

Description	Mouse monoclonal antibody against CD130 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Sf21 derived recombinant human CD130 (Leu24-Glu619)
Target	CD130
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P40189
Expiration Date	12 months from date of receipt.