

Product Datasheet

CCL17 antibody (Biotin) (orb867183)

Description	Mouse monoclonal antibody against CCL17 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	IA
Immunogen	E. coli derived recombinant human full length CCL17 (Ala24-Ser94)
Target	CCL17
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q92583
Expiration Date	12 months from date of receipt.