

Product Datasheet

CCL13 antibody (Biotin) (orb867188)

Description	Mouse monoclonal antibody against CCL13 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	E. coli derived recombinant human full length CCL13 (Phe17-Thr98)
Target	CCL13
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q99616
Expiration Date	12 months from date of receipt.