
Product Datasheet

Cytokeratin 10+13 antibody (orb867241)

Description

Mouse monoclonal antibody to Cytokeratin

Species/Host

Mouse

Reactivity

Canine, Feline, Human, Zebrafish

Conjugation

Unconjugated

Tested Applications

FC, ICC, IHC-Fr, IHC-P, WB

Preservatives

Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

Form/Appearance

Each vial contains 100 ul 1 mg/ml purified monoclonal antibody in PBS containing 0.09% sodium azide.

Storage

Storage: The antibody is shipped at ambient temperature and may be stored at +4°C. For prolonged storage prepare appropriate aliquots and store at or below -20°C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance or the concentration of the product

Note

For research use only

Isotype

IgG2a

Clonality

Monoclonal

Source

DE-K13 is a Mouse monoclonal IgG2a, κ antibody derived by fusion of SP2/0 Mouse myeloma cells with spleen cells from a (BALB/c x B6)F1 Mouse immunized with a cytoskeletal preparation extracted from Human ectocervical epithelium.

Uniprot ID

P13645, P13646

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not