
Product Datasheet

Human Tau-441 / 2N4R (241-380) Protein (orb867321)

Description

Human Tau-441 (241-380), His Tag is

AA Ser 241 - Glu 380 (...)

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Tau

Preservatives50 mM Tris, 150 mM NaCl, 1 mM TCEP,
1 mM EDTA, pH7.5**Form/Appearance**

Powder

Storage

-20°C

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 16.9 kDa. The protein migrates as 22 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence

NP_005901.2

Purity

95%

Source

Human Tau-441 (241-380), His Tag (orb867321) is expressed from E. coli cells. It contains AA Ser 241 - Glu 380 (Accession # P10636-8).

MW

16.9 kDa

Uniprot ID[P10636](#)**NCBI**[NP_005901.2](#)**Expiration Date**

6 months from date of receipt.