
Product Datasheet

Human Latent TGF-beta 3 / Latent TGFB3 Protein (orb867322)

Description

Human Latent TGF-beta 3, His Tag is

(HEK293). It contains AA Leu 24 -...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TGF-beta 3

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 46.8 kDa. The protein migrates as 44-47 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_003230.1

Purity

90%

Source

Human Latent TGF-beta 3, His Tag (orb867322) is expressed from human 293 cells (HEK293). It contains AA Leu 24 - Ser 412 (Accession # P10600-1).

MW

46.8 kDa

Uniprot ID[P10600](#)**NCBI**[NP_003230.1](#)**Expiration Date**

6 months from date of receipt.