

Product Datasheet

Biotinylated Mouse TSLP Protein (orb867327)

Description	Biotinylated Mouse TSLP, His, is expressed from human 293 cells (HEK293). It contains AA Tyr 20 - Gl...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Biotin
Target	TSLP
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His & C-Avi
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 17.6 kDa. The protein migrates as 25-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Tyr 20 - Glu 140
Purity	90%
Source	Biotinylated Mouse TSLP, His, (orb867327) is expressed from human 293 cells (HEK293). It contains AA Tyr 20 - Glu 140 (Accession # Q9JIE6-1).
MW	17.6 kDa
Uniprot ID	Q9JIE6
Expiration Date	6 months from date of receipt.