

Product Datasheet

Mouse TPBG / 5T4 Protein (orb867330)

Description	Mouse TPBG, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 32 - Ser 361 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TPBG
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.7 kDa. The protein migrates as 55-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_035757.2
Purity	95%
Source	Mouse TPBG, His Tag (orb867330) is expressed from human 293 cells (HEK293). It contains AA Ser 32 - Ser 361 (Accession # Q9Z0L0-1).
MW	37.7 kDa
Uniprot ID	Q9Z0L0
NCBI	NP_035757.2
Expiration Date	6 months from date of receipt.