

Product Datasheet

Mouse MERTK / Mer Protein (orb867338)

Description	Mouse MERTK, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gly 19 - Met 497 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	MERTK
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 78.9 kDa. The protein migrates as 100-150 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_032613.1
Purity	95%
Source	Mouse MERTK, Fc Tag (orb867338) is expressed from human 293 cells (HEK293). It contains AA Gly 19 - Met 497 (Accession # Q60805-1).
MW	78.9 kDa
Uniprot ID	Q60805
NCBI	NP_032613.1
Expiration Date	6 months from date of receipt.