

## Product Datasheet

### SCRIPT RT-qPCR SybrMaster UNG (orb867347)

<b>Description</b>	SCRIPT RT-qPCR SybrMaster UNG is designed for quantitative real-time analyses of RNA templates using...
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	2x conc
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Spectroscopic Propertie:</b> $\lambda_{exc}$ 494 nm (bound to DNA), $\lambda_{em}$ 521 nm, (bound to DNA).
<b>Expiration Date</b>	12 months from date of receipt.