

Product Datasheet

Mouse SIGLEC-15 protein (orb867405)

Description	Recombinant Mouse SIGLEC-15 (C-hFc) is produced by Human Cells. The target gene encoding Arg24-Thr26...
Endotoxins	Less than 1 EU/ug
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Human Fc
Note	For research use only
Uniprot ID	Q6ZMC9
Expiration Date	12 months from date of receipt.