

## Product Datasheet

### Human CCL20 protein (orb867425)

<b>Description</b>	Recombinant Human CCL20 (C-hFc) is produced by Human Cells. The target gene encoding Ala27-Met96 is ...
<b>Endotoxins</b>	Less than 1 EU/ug
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-Human Fc
<b>Note</b>	For research use only
<b>Uniprot ID</b>	AAH20698.1
<b>Expiration Date</b>	12 months from date of receipt.