

Product Datasheet

Cynomolgus TREM1 protein (orb867436)

Description	Recombinant Cynomolgus TREM1 (C-hFc) is produced by Human Cells. The target gene encoding Thr21-Pro2...
Endotoxins	Less than 1 EU/ug
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Human Fc
Note	For research use only
Uniprot ID	A0A1D5RC96
Expiration Date	12 months from date of receipt.