

Product Datasheet

Human B7H3-4Ig protein (orb867452)

Description	Recombinant Human B7H3-4Ig/B7-H3b (C-His) is produced by Human Cells. The target gene encoding Gly27...
Endotoxins	Less than 1 EU/ug
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Uniprot ID	Q5ZPR3
Expiration Date	12 months from date of receipt.