

## Product Datasheet

### Mouse CD276 protein (orb867455)

<b>Description</b>	Recombinant Mouse CD276/B7-H3 (C-mFc) is produced by Human Cells. The target gene encoding Glu30-Ala...
<b>Endotoxins</b>	Less than 1 EU/ug
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-Mouse Fc
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q8VE98</b>
<b>Expiration Date</b>	12 months from date of receipt.