
Product Datasheet

cdc25A (phospho-Ser76) antibody (FITC) (orb8677)

Description

Rabbit polyclonal to cdc25A (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human cdc25A around the phosphorylation site of Ser76 [MG(p-S)E]
Target	CDC25A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59kDa
Uniprot ID	P30304
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.