

## Product Datasheet

### Goat Anti-Rat IgG H&L antibody (orb868164)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Goat Anti-Rat IgG H&L antibody   |
| <b>Species/Host</b>    | Goat   |
| <b>Reactivity</b>      | Rat  |
| <b>Conjugation</b>     | Unconjugated   |
| <b>Immunogen</b>       | Native Rat IgG   |
| <b>Preservatives</b>   | 0.01M PBS (pH7.4)  |
| <b>Form/Appearance</b> | Lyophilized powder   |
| <b>Storage</b>         | Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | IgG  |
| <b>Clonality</b>       | Polyclonal   |
| <b>Expiration Date</b> | 12 months from date of receipt.  |