

Product Datasheet

Rabbit Anti-Goat IgG H&L antibody (orb868168)

Description	Rabbit Anti-Goat IgG H&L antibody
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Unconjugated
Immunogen	Native Goat IgG
Preservatives	0.01M PBS (pH7.4)
Form/Appearance	Lyophilized powder
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.