

## Product Datasheet

### Goat Anti-Rabbit IgG H&L antibody (orb868170)

<b>Description</b>	Goat Anti-Rabbit IgG H&L antibody
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Native Rabbit IgG
<b>Preservatives</b>	0.01M PBS (pH7.4)
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.