

Product Datasheet

Mouse Anti-Rabbit IgM/Bio antibody (orb868216)

Description	Mouse Anti-Rabbit IgM/Bio antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	Native Rabbit IgM
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	2.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,IF=1:100-1000,ELISA=1:1000-10000
Expiration Date	12 months from date of receipt.