

## Product Datasheet

# Goat Anti-Mouse IgM antibody (Biotin) (orb868218)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Goat Anti-Mouse IgM antibody (Biotin)   |
| <b>Species/Host</b>        | Goat  |
| <b>Reactivity</b>          | Mouse   |
| <b>Conjugation</b>         | Biotin  |
| <b>Tested Applications</b> | ELISA, IF, IHC-Fr, IHC-P, WB  |
| <b>Immunogen</b>           | Native Mouse IgM  |
| <b>Preservatives</b>       | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.           |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 2.0 mg/ml   |
| <b>Storage</b>             | Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                   |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Dilution Range</b>      | WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,IF=1:100-1000,ELISA=1:1000-10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |