

Product Datasheet

Rabbit Anti-Human IgG F(ab')₂ antibody (Biotin) (orb868219)

Description	Rabbit Anti-Human IgG F(ab') ₂ antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	Human IgG F(ab') ₂ fragment
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,IF=1:100-1000,ELISA=1:1000-10000
Expiration Date	12 months from date of receipt.