

## Product Datasheet

# Rabbit Anti-Rat IgG H&L antibody (Biotin) (orb868237)

<b>Description</b>	Rabbit Anti-Rat IgG H&L antibody (Biotin)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IF, IHC-Fr, IHC-P, WB
<b>Immunogen</b>	Native Rat IgG
<b>Preservatives</b>	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	2.0 mg/ml
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB=1:1000-10000,IHC-P=1:100-500,IHC-F=1:100-1000,IF=1:100-1000,ELISA=1:1000-10000
<b>Expiration Date</b>	12 months from date of receipt.