

Product Datasheet

Donkey Anti-Goat IgG H&L antibody (PE-Cy3) (orb868510)

| | |
|----------------------------|---|
| Description | Donkey Anti-Goat IgG H&L antibody (PE-Cy3) |
| Species/Host | Donkey |
| Reactivity | Goat |
| Conjugation | PE/Cy3 |
| Tested Applications | FC |
| Immunogen | Native Goat IgG |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 2.0 mg/ml |
| Storage | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | IF=1:100-1000,Flow-Cyt=1:100-1000 |
| Expiration Date | 12 months from date of receipt. |