

Product Datasheet

Mouse Anti-Rat IgG H&L antibody (PE-Cy7) (orb868529)

Description	Mouse Anti-Rat IgG H&L antibody (PE-Cy7)
Species/Host	Mouse
Reactivity	Rat
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	Native Rat IgG
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	2.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF=1:100-1000,Flow-Cyt=1:100-1000
Expiration Date	12 months from date of receipt.