

Product Datasheet

IRDye 800CW Goat anti-Rabbit IgG (H+L) antibody (orb868544)

| | |
|----------------------------|---|
| Description | IRDye 800CW Goat anti-Rabbit IgG (H+L) antibody |
| Species/Host | Goat |
| Reactivity | Rabbit |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Native Rabbit IgG |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.01% Sodium Azide and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1.0 mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Polyclonal |
| Dilution Range | WB=1:10000-30000 |
| Expiration Date | 12 months from date of receipt. |