

Product Datasheet

Mouse Anti rabbit IgG mAb, HRP conjugated antibody (orb868551)

Description	Mouse Anti rabbit IgG mAb, HRP conjugated antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Native Rabbit IgG
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB=1:1000-10000,ELISA=1:1000-10000
Expiration Date	12 months from date of receipt.