

## Product Datasheet

# Rabbit Anti-Goat IgG H&L antibody (PE-Cy3) (orb868579)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit Anti-Goat IgG H&L antibody (PE-Cy3)                              |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Goat  |
| <b>Conjugation</b>         | PE/Cy3  |
| <b>Tested Applications</b> | FC  |
| <b>Immunogen</b>           | Native Goat IgG   |
| <b>Preservatives</b>       | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 2.0 mg/ml   |
| <b>Storage</b>             | Store at -20°C for one year. Avoid repeated freeze/thaw cycles.         |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Dilution Range</b>      | IF=1:100-1000,Flow-Cyt=1:100-1000                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |