

## Product Datasheet

### Goat Anti-rabbit IgG H&L/PE-Cy7 antibody (orb868598)

<b>Description</b>	Goat Anti-rabbit IgG H&L/PE-Cy7 antibody
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	PE/Cy7
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Native rabbit IgG
<b>Preservatives</b>	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	2.0 mg/ml
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified by Protein G, nonspecific adsorbed
<b>Dilution Range</b>	IF=1:100-1000,Flow-Cyt=1:100-1000
<b>Expiration Date</b>	12 months from date of receipt.