

Product Datasheet

Rabbit Anti-Bovine IgG H&L antibody (APC) (orb868705)

Description	Rabbit Anti-Bovine IgG H&L antibody (APC)
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	APC
Tested Applications	FC
Immunogen	Native Bovine IgG
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	2.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF=1:100-1000,Flow-Cyt=1:100-1000
Expiration Date	12 months from date of receipt.