

Product Datasheet

Rabbit Anti-Human IgG F(ab')₂ antibody (PE-Cy5) (orb868774)

Description	Rabbit Anti-Human IgG F(ab') ₂ antibody (PE-Cy5)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5
Tested Applications	FC
Immunogen	Human IgG F(ab') ₂ fragment
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF=1:100-1000,Flow-Cyt=1:100-1000
Expiration Date	12 months from date of receipt.