

Product Datasheet

Rabbit Anti-Human IgG F(ab')₂ antibody (PE-Cy7) (orb868777)

| | |
|----------------------------|---|
| Description | Rabbit Anti-Human IgG F(ab') ₂ antibody (PE-Cy7) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | PE/Cy7 |
| Tested Applications | FC |
| Immunogen | Human IgG F(ab') ₂ fragment |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1.0 mg/ml |
| Storage | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | IF=1:100-1000,Flow-Cyt=1:100-1000 |
| Expiration Date | 12 months from date of receipt. |